Medovations SafeGuide™ Over the Guidewire Esophageal Dilatation System

*Directions for use and cleaning.
Read carefully prior to use.*

Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

102 E. Keefe Ave.
Milwaukee, WI 53212 USA
800.558.6408
414.265.7620
Fax 414.265.7628
medo@medovations.com
www.medovations.com

Authorized Representative in EU:
MEDIMARK®EUROPE
BP2332. 38033 GRENOBLE CEDEX
FRANCE

Made in the U.S.A.

Medovations, Inc.
102 E. Keefe Ave.
Milwaukee, WI 53212 USA
800.558.6408
414.265.7620
Fax 414.265.7628
medo@medovations.com
www.medovations.com

MEDOVARIONS SAFEGUIDE™ OVER THE GUIDEWIRE DILATATION SYSTEM

Directions for Use and Cleaning.

Read carefully prior to use.

DESCRIPTION

The SAFEGUIDE™ OVER THE GUIDEWIRE ESOPHAGEAL DILATATION SYSTEM includes:

1. A set of flexible, tapered, polyvinylchloride (PVC) dilators that have a central, longitudinal lumen for over-the-guidewire placement.
2. A marked, stainless steel Spring Tip Guidewire, 210cm long with a removable end cap.
3. Cleaning accessories such as a brush and disposable cleaning adaptors.
4. A storage case.

ALL OF THE ABOVE COMPONENTS MAY BE PURCHASED SEPARATELY.

INTENDED USE

For dilatation of upper esophageal webs, lower esophageal rings, caustic strictures, peptic esophageal strictures, and temporary ease of esophageal carcinoma.

The SafeGuide™ Over the Guidewire Dilators are designed to be used over a pre-positioned Medovations SafeGuide™ Guidewire, product number 1214-02.

INSPECTION – PRIOR TO USE

Inspect new dilators and accessories for damage from shipping and handling. Inspect older dilators for wear. Worn dilators should NOT be used. Signs of wear can include: enlarged distal tip lumen due to repeated abrasions over the guidewire, surface cracks and crazing, cuts, and extremely faded print.

There is not an expiration date for the dilators because fitness for use is not age dependent. Fitness for use is determined by inspection for wear.

Visually inspect the guidewire for signs of wear or degradation such as kinks, bends, or breaks. Worn guidewires should NOT be used as they may not perform as expected.

CONTRAINDICATIONS

Contraindications include those specific to upper GI endoscopy.

Contraindications to dilation include, but are not limited to: uncooperative patient; asymptomatic strictures; inability to advance the dilator through the strictured area; coagulopathy; known or suspected perforation; severe inflammation or scarring near the dilation site, recent myocardial infarction, active ulcer and severe cervical arthritis.

POTENTIAL COMPLICATIONS

Potential complications associated with upper gastrointestinal endoscopy and esophageal dilation include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

INSTRUCTIONS FOR USE

1. Perform screening endoscopy to identify strictured area.
2. Remove the protective end cap from the guidewire. Introduce the guidewire, spring tip end first, through the accessory channel of endoscope. A lubricant will help the guidewire glide through smoothly. Advance the guidewire until the spring tip is endoscopically visualized well beyond tip of scope and beyond the strictured area. The guidewire can also be externally monitored using the marking bands. The guidewire is marked in 20cm increments to help determine the location of the spring tip. The markings are groups of bands, beginning distally with a group of two (2) bands 40 cm from the tip. Each band represents 20 cm of distance to the spring tip, ending in a grouping of seven (7) bands 140 cm from the tip. A simple formula to measure the distance from a specific group of bands to the spring tip is:
$$(\text{number of bands}) \times 20 \text{ cm} = \text{distance in centimeters to the spring tip}$$
Measurements are taken from the most distal band in each group
CAUTION: When placing the guidewire, discontinue advancement of the wire if resistance is met.
3. When the guidewire is in position well beyond strictured area, slowly begin to withdraw endoscope in 5-10cm increments while simultaneously advancing guidewire in 5-10cm increments to ensure guidewire remains in position. CAUTION: Continuous fluoroscopic monitoring of the guidewire is essential in order to ensure it remains in proper position.
4. When endoscope is removed completely, fluoroscopically confirm that guidewire has not been displaced.
Note: The guidewire position can be externally monitored using the marking bands (as described in 2 above) in reference to the dental arch, both before and during dilatation.

Medovations, Inc.
102 E. Keefe Ave.
Milwaukee, WI 53212 USA
800.558.6408
414.265.7620
Fax 414.265.7628
medo@medovations.com
www.medovations.com

5. Generously lubricate dilator and advance it over pre-positioned guidewire to strictured area.
Note: Medovations SafeGuide™ Dilators have two (2) series of numbers printed on each dilator on opposing sides. One series measures the distance to the distal tip, identified as FROM TIP (American). The other series measures the distance to the stricture, identified as FROM MAX O.D. (European). These printed distance markings can be used as a guide to avoid over-extending into the gastric cavity.
6. Proceed with esophageal dilation. Determine the SafeGuide™ dilator to be used by the size of the strictures, starting with a smaller French sized dilator and gradually increasing the French size. Hold the dilator at the blunt, proximal end and place the lumen at the distal tapered tip over the end of the guidewire. Slide the dilator along the guidewire while continuing to monitor the guidewire location. Carefully remove dilator after completed dilation.
Note: Successive dilations can be accomplished by repeating steps 5-6.
7. Upon completion of esophageal dilation, remove dilator and guidewire from patient.
8. Clean and prepare dilator for future use following cleaning instructions below.

CLEANING INSTRUCTIONS

CAUTION:

1. During cleaning, inspect integrity and function of device to determine advisability of reuse. If kinks, excessive wear, bends or breaks exist, do not use. Dispose of device per institutional guidelines for biohazardous medical waste.
2. DO NOT USE Betadine to clean the dilators as it may cause discoloration of the material.
3. SafeGuide™ Dilators CANNOT be steam autoclaved. The high temperatures of steam autoclaving will destroy the dilator.
4. DO NOT USE alcohol for cleaning or disinfection.
5. Appropriate gloves, gowns, and other protective wear should be worn when cleaning and sterilizing soiled or used dilators.
6. It is the responsibility of the user to qualify any deviations from the recommended methods of cleaning.

Medovations™ SafeGuide™ Dilators must be thoroughly cleaned prior to use.

PRE-CLEANING (SafeGuide™ Dilators and Spring Tip Guidewires):

1. Immediately after use, soak dilator/ guidewire in a solution of warm water and enzymatic cleaning solution. Note: Refer to the enzymatic cleaning solution instructions for recommended duration of cleaning.
2. Remove gross contaminates on the external surfaces of the dilator/guidewire by hand, using a soft bristle brush, sponge, or wipe, with warm water and enzymatic solution.
3. Clean the inner lumen of dilator using a solution of warm water and enzymatic solution and the Medovations SafeGuide™ Cleaning Brush. The WIRE END of the Cleaning Brush must be inserted first into the lumen of the dilator. Pull the brush through the dilator, brushing the inner lumen of the dilator.
CAUTION: DO NOT PUSH THE CLEANING BRUSH through the dilator with the brush end first. Pushing the Cleaning Brush through the dilator, brush end first, will cause the wire to kink and/or the brush to break.
4. Use the SafeGuide™ Cleaning Adaptors to flush the inner lumen with a solution of warm water and enzymatic cleaner.
5. Thoroughly rinse the dilator/guidewire with clean water. Use the SafeGuide™ Cleaning Adaptors to flush the inner lumen of the dilator with clean water.
6. Pat dry with clean, non-contaminating wipes. Use medical grade (oil-free) forced air to dry the lumen of the dilator.
7. Proceed with high-level disinfection / sterilization of dilator/guidewire.

The SafeGuide™ Guidewire (without the protective end cap) and SafeGuide™ Cleaning Brush may be exposed to steam autoclave, using a standard steam sterilization cycle of 132°C for 15 minutes. The guidewires and cleaning brushes should be individually wrapped for the steam sterilization cycle. After steam sterilization, allow the guidewires/cleaning brushes to completely cool to room temperature (25°C). Do not immerse in cold liquids. Dry thoroughly.

Medovations, Inc.
102 E. Keefe Ave.
Milwaukee, WI 53212 USA
800.558.6408
414.265.7620
Fax 414.265.7628
medo@medovations.com
www.medovations.com

The SafeGuide™ Guidewire (without the protective end cap) and SafeGuide™ Cleaning Brush may also be exposed to a glutaraldehyde soak.

HIGH-LEVEL DISINFECTION / STERILIZATION:

The following disinfection solutions can be used with the Medovations SafeGuide™ Over the Guidewire Esophageal Dilatation System:

- Glutaraldehyde
- Ortho-phthalaldehyde (OPA)
- Peracetic acid
- Hydrogen peroxide

Please refer to the manufacturer's instructions for complete information and precautions for disinfection solutions with any of the above active ingredients.

The dilators have additionally been validated with the following specific brands:

- Cidex® OPA High Level Disinfectant
- Steris Resert™ XL HLD High Level Disinfectant
- Rapicide PA High Level Disinfectant
- Reliance™ EPS Endoscope Processing System
- Medivators Advantage Plus System
- Medivators DSD series of automated endoscope reprocessors
- Steris SYSTEM 1E® Liquid Chemical Sterilant Processing System

Please refer to the manufacturer's instructions for complete information and precautions.

Note: Proper rinsing with clean water and complete drying of the dilator is necessary with any disinfection solution.

1. Completely immerse dilator in a disinfectant. Note: Refer to disinfectant solution instructions for duration of immersion required for high level disinfection or sterilization. If difficulty is encountered with solution penetration into the guidewire lumen, attach suction to the tapered end of dilator or utilize a SafeGuide™ Cleaning Adaptor.
2. Rinse dilator thoroughly with clean, running water. Rinse the guidewire lumen by attaching suction to tapered end of dilator or using a SafeGuide™ Cleaning Adaptor to suction water through this inner guidewire lumen. Note: Proper rinsing with clean water and complete drying of dilator and guidewire lumen is necessary following use of these solutions.

3. Dry external portions for dilator completely with a lint-free cloth. Use medical grade (oil-free) forced air to dry guidewire lumen.

STORAGE

The SafeGuide™ dilators should be stored flat, out of direct light, and at room temperature. The SafeGuide™ Storage Case provides a protected environment.

Medovations, Inc.
102 E. Keefe Ave.
Milwaukee, WI 53212 USA
800.558.6408
414.265.7620
Fax 414.265.7628
medo@medovations.com
www.medovations.com

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